

DSC CONSUMABLES

incorporated

We specialize in the most frequently used DSC sample pans and thoroughly clean every order prior to packaging & shipping. This disciplined process ensures our customers receive contaminant free products upon delivery saving you time and resources



Qty/Pkg: 200 Code: DSC84001

Standard aluminum sample pans compare to TA Instruments #900786.901



Qty/Pkg: 200 Code: DSC84002

Standard aluminum lids compare to TA Instruments #900779.901.



Qty/Pkg: 200 Code: DSC84003

Standard aluminum sample pan/lid set compare to TA Instruments #'s 900786.901/900779.901.



Qty/Pkg: 200 Code: DSC84004

Hermetic aluminum sample pans compare to TA Instruments #900793.901.



Qty/Pkg: 200 Code: DSC84005

Hermetic aluminum lids compare to TA Instruments #900794.901.



Qty/Pkg: 200 Code: DSC84006

Hermetic aluminum sample pan/lid set compare to TA Instruments #'s 900793.901/900794.901.



Qty/Pkg: 100 Code: DSC22001

Aluminum crucible standard, 40ul, w/o pin compare to Mettler.



Qty/Pkg: 100 Code: DSC22002

Aluminum lids compare to Mettler ME 51119871.



Qty/Pkg: 100 Code: DSC22003

Aluminum crucible with lid set compare to Mettler (flat) set ME-00026763.



Qty/Pkg: 100 Code: DSC22004

Aluminum crucible standard, 40ul, w/ pin compare to Mettler.



Qty/Pkg: 100 Code: DSC22005

Aluminum lids compare to Mettler ME 51119871.



Qty/Pkg: 100 Code: DSC22006

Aluminum crucible with lid set compare to Mettler w/pin set ME-00027331.



Qty/Pkg: 100 Code: DSC22013

Aluminum crucible with lid set 100uL, compare to Mettler.



Qty/Pkg: 10 Code: DSC22008

Premium alumina crucible/lid set (70uL) compare to Mettler 00024123.



Qty/Pkg: 10 Code: DSC22009

Premium alumina crucible/lid set (150uL) compare to Mettler 00024124.



Qty/Pkg: 2 Code: DSC22010

Sapphire 6.0mm x 2.6mm crucible/lid set compare to Mettler 30uL.



Qty/Pkg: 2 Code: DSC22011

Sapphire 6.0mm x 4.5mm crucible/lid set compare to Mettler 70uL.



Qty/Pkg: 2 Code: DSC22012

Sapphire 8.0mm x 4.5mm crucible/lid set compare to Mettler 150uL.



Qty/Pkg: 1 Code: DSC72001

Premium alumina crucible compare to Netzsch #399972



Qty/Pkg: 1 Code: DSC72002

Premium alumina lid compare to Netzsch #399973



Qty/Pkg: 1 Code: DSC72003

Premium alumina crucible/lid set compare to Netzsch #399972/399973



Qty/Pkg: 100 Code: DSC72004

Standard aluminum sample pans compare to Netzsch 40uL.



Qty/Pkg: 100 Code: DSC72005

Standard aluminum sample lids compare to Netzsch.



Qty/Pkg: 100 Code: DSC72006

Standard aluminum sample pan/lid set compare to Netzsch 40uL.



Qty/Pkg: 2 Code: DSC72007

Sapphire 6.8mm x 4.0mm crucible/lid set compare to Netzsch 85uL.



Qty/Pkg: 200 Code: DSC64001

Aluminum sample pans compare to Perkin Elmer standard pans.



Qty/Pkg: 200 Code: DSC64002

Aluminum lids compare to Perkin Elmer standard lids.



Qty/Pkg: 200 Code: DSC64003

Standard aluminum sample pan/lid set compare to Perkin Elmer #02190041.



Qty/Pkg: 200 Code: DSC44001

Aluminum sample pans compare to Seiko standard pans.



Qty/Pkg: 200 Code: DSC44002

Aluminum lids compare to Seiko standard lids.



Qty/Pkg: 200 Code: DSC44003

Standard aluminum sample pan/lid set compare to Seiko standard pan/lid set.



Qty/Pkg: 1 Code: DSC44004

5 mm x 5 mm platinum pan.



Qty/Pkg: 1 Code: DSC44005

5 mm x 5 mm alumina crucible w/cover.

Sapphire Crucibles



Qty/Pkg: 2 Code: DSC22010

Sapphire 6.0mm x 2.6mm crucible/lid set compare to Mettler 30uL.



Qty/Pkg: 2 Code: DSC22011

Sapphire 6.0mm x 4.5mm crucible/lid set compare to Mettler 70uL.



Qty/Pkg: 2 Code: DSC22012

Sapphire 8.0mm x 4.5mm crucible/lid set compare to Mettler 150uL.



Qty/Pkg: 2 Code: DSC72017

Sapphire 6.8mm x 4.0mm crucible/lid set compare to Netzsch 85uL.



Qty/Pkg: 2 Code: DSC84818

Sapphire 6.5mm x 4.3mm crucible/lid set compare to TA Instruments 90uL.